

50 ml tubes

CAT#: TCC-TU001-50

QTY: 1 x 25 tubes

STORAGE CONDITIONS: Room temperature

SHIPPING: Room temperature

DESCRIPTION:

These tubes are equipped with caps that have a crosscut permeable filter, preventing airborne particulates from contaminating samples during run cycles.

INSTRUCTIONS OF USE:

These 50ml tubes are used exclusively in the outlet (collection) ports to recover sorted samples. If a sample run is expected to take more than an hour, use a 50 mL tube for collection.

UNITED STATES

1100 Island Drive, Suite 203, Redwood City, CA

JAPAN

7-3-1 Hongo, Bunkyo, Tokyo

THINKCYTE.COM

Copyright © 2023 ThinkCyte, K.K.. ThinkCyte, Ghost Cytometry, and VisionSort are trademarks of ThinkCyte, K.K.. All rights reserved.
For Research Use Only.