

VS Clean

rev.9/24

CAT#:	TCR-CN001-250, TCR-CN001-1000
AMOUNT:	250 mL, 1000 mL
FORM:	Liquid
STERILITY:	0.22 µm sterile filtered under sterile conditions
SHELF LIFE:	Stable for 12 months under proper storage and handling
STORAGE CONDITIONS:	Room Temperature, Light sensitive
SHIPPING:	Room temperature

DESCRIPTION:

VS Clean is an alkaline-based cleansing solution that helps remove nucleic acids and cellular debris from the instrument. It is also used to remove clogs during routine instrument cleaning cycles.

INSTRUCTIONS FOR USE:

Apply **5 mL** of VS Clean solution for each of the two sample inlets of VisionSort for cleaning at the end of each day.

Note: Avoid running VS Clean on a VS Sort Cartridge